



Audiometer ear 3.0

Modern Applications

mobile or integrated

AURITEC 



Audiometer ear 3.0

PC-driven, state-of-the-art Audiometer

- Easy installation due to an integrated soundcard, which is connected to the Notebook or the PC via USB, and easy, automated installation program
- Speech tests which are available on CD, can be easily stored on the hard drive and be accessed from there
- Improved mobility through an optimised use with Notebook, Akkupak and transportation suitcase
- To be used with Windows xp professional and Windows 7



- Integrated into the NOAH-platform; results can also be stored within the device
- Applicable with Laptop or Desktop-PC, mobil or used in the cabin, different loudspeakers available, upgradeable to a Regiometry
- The ear 3.0 can be combined with hearing-aid fitting systems – becoming a complete package for fitting
- Modular design – all options can be upgraded later, for example high frequency audiometry

The logo for AURITEC, featuring the word "AURITEC" in a bold, blue, sans-serif font. To the right of the text is a stylized orange arc that curves around the top right corner of the letters. The background of the entire top section is a light blue grid with faint white lines and a larger, faint white wave-like pattern on the right side.

AURITEC



AURITEC Medizindiagnostische Systeme GmbH

Dernauer Straße 12
D-22047 Hamburg

Telefon: +49 (0)40-69 65 40-0
Telefax: +49 (0)40-69 65 40-10

E-Mail: ear3.0@auritec.de
Internet: www.auritec.de