


# Picture Presentation – VRA Viewer

- › For conditioning of children with pictures on flat screens.
- › Simple handling in connection with PTA.
- › Compatible with every PC driven audiometer, especially AT2000-PED
- › Easy to use, integrated on the audiometry screen
- › Possibility to complete the picture and video database individually.



# Picture Presentation – VRA Viewer

This is a special software for the conditioning of children during audiometrical diagnosis.

The aim is to keep the children's attention for as long as possible in order to gain more time for the different parts of your audiological work.

Children nowadays are used to colorful pictures and videos. That is the reason why AURITEC uses these medias to catch the child's eyes and brain.

The range of pictures has been selected with doctors from German University clinics and also together with experienced audiologists.

You can present the pictures or videos via the cursor keys of the laptop or PC or via a IR remote control.

The format of the presented pictures matches the generally used picture formats (bmp, jpg).

Please also ask for the Looky Box for reinforcement of younger children



## Delivered devices

- Grafic card
- 2 TFT flat screens
- Software VRA 4.0 with pictures and videos

## Option

- Installation

